



Product Datasheet

EGFL7 antibody (orb13387)

Description	EGFL7 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse EGFL7 (31-130/275 aa)
Target	EGFL7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	Q9QXT5
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human ovarian carcinoma tissue ...