


Product Datasheet

Cytokeratin 18 (KRT18) antibody (orb1338701)

Description	Mouse monoclonal antibody to Cytokeratin 18 (KRT18)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 69-372 of human CK18 (NP_000215) produced in E.coli.
Target	KRT18
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	47.9 kDa
Uniprot ID	P05783
NCBI	NM_199187
Dilution Range	WB 1:1000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...