


Product Datasheet

c-Myc (MYC) antibody (orb1338710)

Description	Mouse monoclonal antibody to c-Myc (MYC)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human Myc (NP_002458) produced in E.coli.
Target	MYC
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	50.4 kDa
Uniprot ID	P01106
NCBI	NM_002467
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...