


Product Datasheet

pan Cytokeratin antibody (orb1338713)

Description	Mouse monoclonal antibody to pan Cytokeratin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human epidermal cytokeratins
Target	pan Cytokeratin
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC 1:1000-1:2500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...