


Product Datasheet

Pan-CDH antibody (orb1338714)

Description	Mouse monoclonal antibody to Pan-CDH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide corresponding to the conserved region of human cadherin family.
Target	Pan-CDH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...