


Product Datasheet

BrdU antibody (orb1338740)

Description	Mouse monoclonal antibody to BrdU
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Bromodeoxyuridine coupled to keyhole limpet hemocyanin (KLH).
Target	BrdU
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC 1:150, IF 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...