


Product Datasheet

Carrier-free (BSA/glycerol-free) tdTomato antibody (orb1338793)

Description	Mouse monoclonal antibody to Carrier-free (BSA/glycerol-free) tdTomato
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A tdTomato tagged recombinant protein expressed in E.coli.
Target	tdTomato
Preservatives	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	54.2 kDa
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...