


Product Datasheet

Carrier-free (BSA/glycerol-free) DsRed-Monomer antibody (orb1338803)

Description	Mouse monoclonal antibody to Carrier-free (BSA/glycerol-free) DsRed-Monomer
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This antibody is raised against a recombinant protein of DsRed-Monomer expressed in HEK293 cells.
Target	DsRed-Monomer
Preservatives	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.4 kDa
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...