

Product Datasheet

Carrier-free (BSA/glycerol-free) mRaspberry/AmCyan1 antibody (orb1338812)

Description Mouse monoclonal antibody to Carrier-free (BSA/glycerol-free) mRaspberry/AmCyan1

Species/Host Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen This antibody is raised against a recombinant protein of mRaspberry expressed in HEK293 cells. It recognizes both mRaspberry and AmCyan1 proteins.

Target mRaspberry

Preservatives Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.


Isotype IgG1

Clonality Monoclonal

MW 25.5 kDa

Dilution Range WB: 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...