


Product Datasheet

Carrier-free (BSA/glycerol-free) AcGFP1/PS-CFP2 antibody (orb1338816)

Description	Mouse monoclonal antibody to Carrier-free (BSA/glycerol-free) AcGFP1/PS-CFP2
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This antibody is raised against a recombinant protein of PS-CFP2 expressed in HEK293 cells. It recognizes both PS-CFP2 and AcGFP1 proteins.
Target	AcGFP1
Preservatives	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	26.7 kDa
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...