



Product Datasheet

EIF4E1 (phospho-Ser209) antibody (orb13390)

Description	EIF4E1 (phospho-Ser209) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human eIF4E around the phosphorylation site of Ser209 (SG(p-S)TT)
Target	phospho-eIF4E (Ser209)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P06730
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=0.2ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Rat brain
using
phosph...



IHC-P
analysis of
Mouse brain
using
phos...