


## Product Datasheet

### Phospholipase C beta 4 (PLCB4) Human Over-expression Lysate (orb1339125)

|                        |                                                                                          |
|------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of phospholipase C, beta 4 (PLCB4), transcript variant 3 |
| <b>Target</b>          | PLCB4                                                                                    |
| <b>Note</b>            | For research use only                                                                    |
| <b>Uniprot ID</b>      | <b>Q15147</b>                                                                            |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                           |



Western blot analysis of Left: Cell lysa...