


## Product Datasheet

### hnRNP Q (SYNCRIP) Human Over-expression Lysate (orb1339797)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of synaptotagmin binding, cytoplasmic RNA interacting protein (SYNCRIP) |
| <b>Target</b>          | SYNCRIP                                                                                                 |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>Q59GL1</b> , <b>O60506</b>                                                                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |



Western blot analysis of Left: Cell lysa...