



Product Datasheet

FOXJ1 antibody (orb13412)


Description	FOXJ1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human FOXJ1 (161-260/421 aa)
Target	HFH4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	Q92949
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500, Flow-Cyt=0.2ug/test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P image of mouse testis tissue using...



Western blot analysis of testis Mouse ly...



Immunohistochemical staining of paraffin...