

## Product Datasheet

# Chk1 (CHEK1) Human Over-expression Lysate (orb1341333)

|                        |                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of<br>CHK1 checkpoint homolog (S.<br>pombe) (CHEK1), transcript variant 3 |
| <b>Target</b>          | CHEK1                                                                                                     |
| <b>Note</b>            | For research use only                                                                                     |
| <b>Uniprot ID</b>      | <b>O14757</b> , <b>B4DT73</b>                                                                             |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                            |