


## Product Datasheet

### Insulin Receptor (INSR) Human Over-expression Lysate (orb1341638)

|                        |                                                                                  |
|------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of insulin receptor (INSR), transcript variant 1 |
| <b>Target</b>          | INSR                                                                             |
| <b>Note</b>            | For research use only                                                            |
| <b>Uniprot ID</b>      | <b>P06213</b>                                                                    |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                   |



Western blot analysis of Left: Cell lysa...