


## Product Datasheet

# HTR2C Human Over-expression Lysate (orb1341923)

|                        |                                                                                        |
|------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of 5-hydroxytryptamine (serotonin) receptor 2C (HTR2C) |
| <b>Target</b>          | HTR2C                                                                                  |
| <b>Note</b>            | For research use only                                                                  |
| <b>Uniprot ID</b>      | <b>P28335</b>                                                                          |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                         |



Western blot analysis of  
Left: Cell lysa...