


## Product Datasheet

### PDE3A Human Over-expression Lysate (orb1341952)

|                        |                                                                                 |
|------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of phosphodiesterase 3A, cGMP-inhibited (PDE3A) |
| <b>Target</b>          | PDE3A                                                                           |
| <b>Note</b>            | For research use only                                                           |
| <b>Uniprot ID</b>      | <a href="#">Q14432</a>                                                          |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                  |



Western blot analysis of Left: Cell lysa...