


## Product Datasheet

### TIRAP Human Over-expression Lysate (orb1343288)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of toll-interleukin 1 receptor (TIR) domain containing adaptor prote... |
| <b>Target</b>          | TIRAP                                                                                                   |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>A0A024R3M4, P58753</b>                                                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |



Western blot analysis of  
Left: Cell lysa...