


## Product Datasheet

### C3orf43 (SMCO1) Human Over-expression Lysate (orb1343680)

|                        |                                                                                 |
|------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of chromosome 3 open reading frame 43 (C3orf43) |
| <b>Target</b>          | SMCO1                                                                           |
| <b>Note</b>            | For research use only                                                           |
| <b>Uniprot ID</b>      | <b>Q147U7</b>                                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                  |



Western blot analysis of Left: Cell lysa...