

## Product Datasheet

### DOPA Decarboxylase (DDC) Human Over-expression Lysate (orb1343814)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of dopa decarboxylase (aromatic L-amino acid decarboxylase) (DDC), t... |
| <b>Target</b>          | DDC                                                                                                     |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <a href="#">A0A0S2Z3N4</a> , <a href="#">Q53Y41</a> , <a href="#">P20711</a>                            |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |

