


## Product Datasheet

# SPAG11A Human Over-expression Lysate (orb1343832)

|                        |                                                                           |
|------------------------|---------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of sperm associated antigen 11A (SPAG11A) |
| <b>Target</b>          | SPAG11A                                                                   |
| <b>Note</b>            | For research use only                                                     |
| <b>Uniprot ID</b>      | <b>A0A0A0MR37, Q6PDA7</b>                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                            |



Western  
blot  
analysis of  
Left: Cell  
lysa...