


## Product Datasheet

### GAGE13 Human Over-expression Lysate (orb1344168)

|                        |                                                          |
|------------------------|----------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of G antigen 13 (GAGE13) |
| <b>Target</b>          | GAGE13                                                   |
| <b>Note</b>            | For research use only                                    |
| <b>Uniprot ID</b>      | <b>Q4V321</b>                                            |
| <b>Expiration Date</b> | 6 months from date of receipt.                           |



Western  
blot  
analysis of  
Left: Cell  
lysa...