


## Product Datasheet

### IRF1 Human Over-expression Lysate (orb1344899)

|                        |                                                                          |
|------------------------|--------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of interferon regulatory factor 1 (IRF1) |
| <b>Target</b>          | IRF1                                                                     |
| <b>Note</b>            | For research use only                                                    |
| <b>Uniprot ID</b>      | <b>Q6FHN8, P10914</b>                                                    |
| <b>Expiration Date</b> | 6 months from date of receipt.                                           |



Western blot analysis of  
Left: Cell lysa...