



Product Datasheet

GST tag antibody (orb13452)

Description	GST tag antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length GST tag protein
Target	GST tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Transformed
(GST-tag or EV)
E. coli cell...



Western
blot: Recombinant
GST Protein; Pri...