


## Product Datasheet

### Transition Protein 1 (TNP1) Human Over-expression Lysate (orb1345421)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of transition protein 1 (during histone to protamine replacement) (T... |
| <b>Target</b>          | TNP1                                                                                                    |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>P09430, Q4ZG82</b>                                                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |



Western blot analysis of Left: Cell lysa...