


Product Datasheet

ERG antibody (orb13456)

Description	ERG antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human HERG (1001-1159/1159 aa)
Target	ERG/KCNH2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	127 kDa
Uniprot ID	Q12809
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Rsc96 cell u...