


## Product Datasheet

### PRMT3 Human Over-expression Lysate (orb1346694)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of protein arginine methyltransferase 3 (PRMT3), transcript variant ... |
| <b>Target</b>          | PRMT3                                                                                                   |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>O60678</b>                                                                                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |



Western  
blot  
analysis of  
Left: Cell  
lysate...