


## Product Datasheet

### SF3A1 Human Over-expression Lysate (orb1346752)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of splicing factor 3a, subunit 1, 120kDa (SF3A1), transcript variant... |
| <b>Target</b>          | SF3A1                                                                                                   |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>A0A024R1K8, Q15459</b>                                                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |



Western blot analysis of  
Left: Cell lysa...