



Product Datasheet

Haptoglobin antibody (orb13477)

Description	Haptoglobin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Haptoglobin beta chain (251-350/406 aa)
Target	Haptoglobin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27/43 kDa
Uniprot ID	P00738
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of human uterine cancer tissue (Ha...)



Immunohistochemical analysis of paraffin...