



Product Datasheet

Syndecan 2 antibody (orb13481)


Description	Syndecan 2 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Syndecan 2. Please contact us for the exact immunogen sequence. The peptide is available as orb12911.
Target	Syndecan 2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	22 kDa
Uniprot ID	P34741 , P34900 , P43407
Entrez	6383
Dilution Range	WB: 1:200-1000, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of guinea pig pancreas tissue



IF image of guinea pig pancreas tissue



IHC-P staining of rat cerebellum tissue