



Product Datasheet

SDHA antibody (orb134901)

Description	Rabbit polyclonal to SDHA. Flavoprotein (FP) subunit of succinate dehydrogenase (SDH) that is involv...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SDHA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72692
Uniprot ID	P31040
NCBI	NP_004159.2
Dilution Range	IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...