



## Product Datasheet

### ECE2 antibody (orb134903)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to ECE2. Converts big endothelin-1 to endothelin-1. Also involved in the processin...                                                                      |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                    |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>Target</b>              | This ECE2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 551-579 amino acids from the Central region of human ECE2.            |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>O60344</b>                                                                                                                                                                |
| <b>NCBI</b>                | <b>NP_001093590.1, NP_115707.2, NP_001032401.1, NP_001093591.1, NP_055508.3</b>                                                                                              |
| <b>Dilution Range</b>      | IHC-P: 1:25, WB: 1:1000                                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Immunohistochemical analysis of formalin...