



Product Datasheet

Calbindin antibody (orb134925)

Description	Rabbit Polyclonal antibody to Calbindin
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CALB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30025
Uniprot ID	P05937
NCBI	NP_004920.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...