


Product Datasheet

SGCG antibody (orb134929)

Description	Rabbit polyclonal to SGCG. Component of the sarcoglycan complex, a subcomplex of the dystrophin-glyc...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SGCG
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32379
Uniprot ID	Q13326
NCBI	NP_000222.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...

