

## Product Datasheet

# CEP55 Human Over-expression Lysate (orb1349401)

|                        |                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of centrosomal protein 55kDa (CEP55), transcript variant 1 |
| <b>Target</b>          | CEP55                                                                                      |
| <b>Note</b>            | For research use only                                                                      |
| <b>Uniprot ID</b>      | <b>Q53EZ4</b>                                                                              |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                             |