

Product Datasheet

Cbwd1 antibody (orb134944)

Description	Rabbit polyclonal to Rat Cbwd1
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Cbwd1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43931
Uniprot ID	Q99MB4
NCBI	NP_598219.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

