

Product Datasheet

Cebpd antibody (orb134946)

Description	Rabbit polyclonal to Rat Cebpd. C/EBP is a DNA-binding protein that recognizes two different motifs:...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Cebpd
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28600
Uniprot ID	Q03484
NCBI	NP_037286.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

