

Product Datasheet

Tcfap2d antibody (orb134947)

Description	Rabbit polyclonal to Mouse Tcfap2d. Sequence-specific DNA-binding protein that interacts with induci...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Tfap2d
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49548
Uniprot ID	Q91ZK0
NCBI	NP_694794.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

