



Product Datasheet

BCL2 antibody (orb134964)

Description	Mouse monoclonal antibody to BCL2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	BCL2L1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	26049
Uniprot ID	Q07817
NCBI	NP_001182.1 , NP_612815.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant pro...