

Product Datasheet

BAIAP2 antibody (orb134984)

Description	Rabbit polyclonal to BAIAP2. Adapter protein that links membrane-bound small G-proteins to cytoplasm...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	BAIAP2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	60868
Uniprot ID	Q9UQB8
NCBI	NP_006331.1 , NP_059345.1 , NP_059344.1 , NP_001138360.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

