

Product Datasheet

STAT1 antibody (orb134986)

Description	Rabbit polyclonal to STAT1. Signal transducer and transcription activator that mediates cellular res...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	STAT1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	87335
Uniprot ID	P42224
NCBI	NP_009330.1 , NP_644671.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

