

Product Datasheet

USP22 antibody (orb134987)

Description	Rabbit polyclonal to USP22. Histone deubiquitinating component of the transcription regulatory histo...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	USP22
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	59961
Uniprot ID	Q9UPT9
NCBI	NP_056091.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

