


Product Datasheet

alpha 2 IL13 Receptor antibody (orb13499)

Description	alpha 2 IL13 Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human IL-13Ra2 (81-180/380 aa)
Target	IL13RA2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	Q14627
Dilution Range	WB=1:200-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Mouse Testis; Lane 2:...