


Product Datasheet

TAS2R13 Human Over-expression Lysate (orb1350678)

Description	Transient overexpression lysate of taste receptor, type 2, member 13 (TAS2R13)
Target	TAS2R13
Storage	Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
MW	35.1 kDa
Uniprot ID	Q9NYV9
Expiration Date	6 months from date of receipt.



Western blot validation of overexpressio...