



Product Datasheet

Inhibin beta B antibody (orb13510)


Description	Inhibin beta B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Inhibin beta B (251-350/407 aa)
Target	Inhibin beta B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13/45 kDa
Uniprot ID	P09529
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Brain (Mouse) L...



Western blot analysis of A431(Human) Cel...



Immunohistochemical analysis of rat ovar...