



Product Datasheet

BTC antibody (orb135143)


Description	Rabbit polyclonal antibody to BTC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BTC(NP_001720.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	P35070
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



IHC-P staining of human lung tissue using the BTC antibody.



Immunohistochemical analysis of formalin-fixed tissue.



Western blot analysis of extracts of various cell lines.