



Product Datasheet

Integrin Alpha V + Beta 1 antibody (orb13515)

Description	Integrin Alpha V + Beta 1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic Mixed peptides 915-948/1048(human Integrin α V)+ 21-70/798(human Integrin β 1)
Target	Integrin Alpha V + Beta 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	114 kDa
Uniprot ID	P06756
Dilution Range	WB=1:500-2000, Flow-Cyt=1 μ g /test, IHC-P=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Sp2/O Lys...



Flow cytometric analysis of Blank contro...