



Product Datasheet

CDK9 antibody (orb135243)


Description	Rabbit polyclonal antibody to CDK9
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CDK9(NP_001252.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40KDa/55KDa
Uniprot ID	P50750
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Western blot analysis of extracts of various cell lines.



Immunofluorescence analysis of U2OS cells.