



Product Datasheet

KLK4 antibody (orb13527)

Description	KLK4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human KLK4 (101-200/254 aa)
Target	KLK4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	Q9Y5K2
Expiration Date	12 months from date of receipt.



IHC-P of human lung cancer using KLK4 an...



Western Blot of mouse brain tissue using...